


High-throughput preparation of specimens for TEM

Michael J. Fasolka, Thomas H. Epps and Kristen Roskov



Polymers
Division

NIST

National Institute of Standards and Technology
Technology Administration, U.S. Department of Commerce



Transmission Electron Microscopy (TEM)

- A primary tool for nanomaterials analysis
- Applicable to a wide range of materials
 - Polymers, metals, ceramics, composites
- Near atomic resolution
- Well understood contrast mechanisms
- Reliable instrumentation




Motivation



Preparation of TEM specimens is difficult and time consuming

Polymeric specimens must be less than ≈ 100 nm thick

- Bulk Specimens
 - Cryogenic ultramicrotomy
 - Slow, one at a time
- Film specimens
 - Sacrificial substrate: float, etch
 - Transparent substrates (SiN)
 - Annealing is a problem
 - No choice of substrate material





A high-throughput method for preparing polymer TEM specimens:


- Thin films
- Use with combinatorial libraries
- Parallel preparation of many specimens
- Flexible
- Amenable to automation

Peel Off Technique



- Originally developed for TEM "detachment replicas"
- Adapted to remove thin films from a substrate

1. Deposit film specimen on substrate
2. Coat film with 10 nm evaporated carbon
3. Pour 25 wt% polyacrylic acid (PAA) solution over specimen
4. Dry PAA layer
5. Peel or "pop" PAA/carbon/specimen from substrate: Razor lift or LN₂ quench
6. Float multilayer on water pool to dissolve PAA
7. Retrieve specimen with TEM grid
8. Voila!




P. H. Geil, in *Polymer Single Crystals* (Wiley Interscience, New York, 1963).
M. J. Fasolka, et al., *Physical Review Letters*, 79, 3018 (1997).


Viability and Advantages of Peel Method




- Films remain intact
- Enables analysis of "buried" film surface
- Enables cross-sectional analysis
- Carbon film adds stability to beam-damaged materials



AFM of bottom surface of BC film peeled from corrugated substrate



↑
Plan and cross-sectional
TEM micrographs of
peeled BC films →




M. J. Fasalca, et al., Physical Review Letters, 79, 3018 (1997).

Strategy for Hi-Throughput Peel Method



1) Fabrication of combinatorial thin film library




Sector Spin Coating*




NCMC Film Methods

Carbon Coat
→

2) Deposit array of PAA droplets




3) Simultaneous removal of droplet/film specimens



Dissolve PAA
→

4) Retrieval of floating film specimens




*de Gans et al, J. COMBI CHEM 7 (6): 952-957 2005

Demonstration of specimen removal



PAA droplets on polystyrene film thickness gradient library on Si





Immersion in pool of LN₂



Droplets "pop" from substrate

PAA/Film specimens harvested from LN₂



Simultaneous specimen removal via LN₂ works...


...but how do we "register" the harvested specimens in library order?

Library registry strategies




1) Pre-labeling of PAA droplets followed by library reconstruction

- Emboss soft PAA droplets with numbers
- Color PAA with gradient of soluble dyes



2) Droplet removal into individual wells

- Spacing of droplets matches well spacing
- Wells filled with water for PAA dissolution




Process	Droplets removed (of 4)
LN ₂	0-1
LN ₂ , sonication	1-2
Sonication, LN ₂ , Sonication	2+

Library Registration Strategies



3) Transfer of droplets to another substrate

- Example: damp glass slide that adheres to PAA




LN₂



Demonstration: 100% film removal with registry



Summary and Future Work



TEM of poly(styrene-b-isoprene)
block copolymer film harvested
from PAA droplet
Thanks to Dr. Brian Berry

Summary:

*Demonstrated a strategy for HT
preparation of TEM specimens*

- Adapted "peel off" method for thin films
 - Samples viable for TEM
 - Amenable to combi libraries
 - Parallel processing
- Evaluated strategies for maintaining registry of library specimens

Future Work:

- Optimize registry methods
- Determine sensitivity to annealing, adhesion, film thickness, specimen type



High-Throughput Facing of Bulk Polymer Specimens


Goal:

- Method for simultaneous *cryogenic fracture facing* of multiple bulk polymer specimens for surface analysis such as AFM and SEM



Experimental Design



- Cryogenic fracture suggests a facing route amenable to multiple sample processing



1. A notched bulk specimen is chilled in LN_2
2. A swift blow to the specimen creates a smooth fracture face amenable to SPM, SEM etc.




AFM of PS PEO Blend - cryofractured

Parallel CryoFracture Technique




By embedding an array of bulk specimens in a notched epoxy matrix, multiple specimen facing could be achieved




Array of 9 different polymer specimens embedded in a notched epoxy capsule

AFM, SEM....



Cryogenic fracture of the epoxy capsule results in facing of the embedded array.



Challenges and Milestones



Challenges to be met:

- Fracture of multiple epoxy-embedded samples not yet demonstrated (single samples have been)
- Multiple sample embedding scheme and capsule design
- Cooling protocol for wide range of specimen T_G s
- Fracture device (Charpy impact tester?)

Deliverables for Focus Project:

- Demonstration of multiple sample fracture
- Embedding protocol and capsule design
- Cooling Protocol
- Demonstration study
- Fracture device/protocol

Impact Areas and Contact Info



- **Impact areas:**
 - Extruded products
 - Polymer Blends
 - Composites and Filled Systems
 - Nanostructured Materials

- **Contact for further discussion:**

Michael J. Fasolka
mfasolka@nist.gov
(301) 975-8526